

OZONOSAN gamma



The latest development from Hänsler Medical for all forms of ozone therapy.

The OZONOSAN gamma is the latest product innovation from Hänsler Medical. It combines proven OZONOSAN quality with remarkable design and functionality.

With an adjustable ozone concentration between 1 and $\leq 80 \mu\text{g/ml}$, all common treatments such as local ozone/oxygen injections, ozonized water (with WOG), major and minor autohemotherapy, rectal insufflation and bag treatment can be performed.

Product highlights



Easy handling via colored touchscreen, supporting multiple language options.

As a tabletop unit, ideal for applications in small and medium-sized medical practices.

Powerful catalyst for safe ozone destruction.

Low maintenance with above-average product service life.

Developed and manufactured in Germany.

Technical specification	OZONOSAN gamma
Device setup	Tabletop unit
Mode of operation	Normobaric infusion technique
Ozone concentration	1 - 80 $\mu\text{g/ml}$
Flow rate	0.9 l/min
Operating pressure	1 bar
Vacuum pump	Up to -0.7 bar
Water treatment with ozone	OZONOSAN WOG optionally available
Dimensions [W x H x D]	42 x 17 x 39 cm
Weight	12 kg
Electrical data	100 - 240 V AC; 50 - 60 Hz; 1.2 A; Protection class I
Classification	IIb medical device according to EU Directive 93/42/EEC
Warranty	2 years; 5 years for ozone tube and catalyst

Hänsler Medical - Only trust the original.

Form of therapy	OZONOSAN gamma
Major autohemotherapy	●
Minor autohemotherapy	●
Rectal insufflation	●
Bag gassing	●
Low pressure ozone gas application	●
Ozone injections	●
Ozonized water	optionally



Ozone therapy treatment sets

OZONOSAN Set 4 Plus for the major autohemotherapy



- 1 OZONOSAN vacuum bottle 250 ml, microbubble system, 12 ml sodium citrate as anticoagulant, latex-free
- 1 Blood transfusion set with steel tip
- 1 Perfusion set optionally with 0.8 or 1.1 mm needle
- 1 Disposable syringe 50 ml, siliconized
- 1 Transfer filter set with cannula and filter (0.2 µm) for ozone introduction

OZONOSAN Set Rectal insufflation with dosing bellow



- 1 Dosing bellow
- 5 Ozone reservoirs with valve
- 50 Special disposable catheters 18 Ch.
- 50 Infusion lines 75 cm
- 50 Bacteria filters

Hänsler Medical GmbH

Nordring 8, 76473 Iffezheim, Germany

+49 (0) 7229 3046 0

@ info@haensler-medical.com

www.haensler-medical.com

