



Hänsler Medical Ozone Generators in comparison.



coming soon



| | OZONOSAN alphaplus | OZONOSAN photonik | OZONOSAN boardcase | OZONOSAN gamma | OZOMED smartline |
|----------------------------|--|---------------------------------------|---------------------------------------|---|--|
| Device setup | Stand-alone, mobile | Tabletop unit | Portable suitcase with wheels | Tabletop unit | Tabletop unit |
| Mode of operation | Normobaric infusion technique | | | | |
| Ozone concentration | 1 - 100 µg/ml | 2 - 75 µg/ml | 1 - 50 µg/ml | 1 - 80 µg/ml | 1 - 80 µg/ml |
| Ozone measurement | Continuous, photometric measurement of the ozone concentration. | | | | |
| Flow rate | 0.5 / 2 / 10 l/min | 0.5 - 10 l/min | 1 - 12 l/min | 0.9 l/min | 0.9 l/min |
| Operating pressure | 0.5 bar | 0.5 bar | 0,5 bar | 1 bar | 1 bar |
| Vacuum pump | Optionally, up to -0.8 bar adjustable | Optionally, up to -0.8 bar adjustable | Optionally, up to -0.8 bar adjustable | Up to -0.7 bar adjustable | Up to -0.7 bar adjustable |
| Ozone extraction | Safety 3-way valve for syringe filling, continous removal and disposal of not used ozone | | | Two separate ozone valves | |
| Water treatment with ozone | 3l reservoir, 1l preperation with ozone | OZONOSAN WOG optionally available | | OZONOSAN WOG optionally available | OZOMED WOG optionally available |
| Dimension [WxHxD] | 45 x 85 x 45 cm | 20 x 56 x 37,5 cm | 55 x 19 x 40 cm | 42 x 17 x 39 cm | 42 x 19 x 39 cm |
| Weight | 38 kg | 20 kg | 17,5 kg incl. oxygen cylinder | 12 kg | 9.8 kg |
| Electrical data | 230 V AC; 50 Hz; 1.2 A; Protection class I | | | 100-240 V AC; 50-60 Hz; 1.2 A; Protection class I | 100-240 V AC; 50-60 Hz; 1.2 A; Protection class I |
| Classification | IIa medical device according to EU Directive 93/42/EEC | | | | IIb medical device according to EU Directive 93/42/EEC |
| Warranty | 2 years; 5 years for ozone tube and catalyst | | | | 2 years in Europe |

Hänsler Medical GmbH - Only trust the original.