

The OZOMED smartline is characterized by its highly compact design and excellent performance. With an adjustable ozone concentration between 1 and $\leq 80 \mu\text{g/ml}$, all common ozone treatments can be performed. The ozonization of water is possible with the optionally available water ozonization unit.

Product highlights

Easy handling via colored touchscreen, supporting multiple language options.

As a tabletop unit, ideal for applications in small and medium-sized medical practices.

Powerful catalyst for safe ozone destruction.

Low maintenance with above-average product service life.

Developed and manufactured in Germany.



| Technical specification | OZOMED smartline |
|----------------------------|--|
| Device setup | Tabletop unit |
| Mode of operation | Normobaric infusion technique |
| Ozone concentration | 1 - 80 $\mu\text{g/ml}$ |
| Flow rate | 0.9 l/min |
| Operating pressure | 1 bar |
| Vacuum pump | Up to -0.7 bar |
| Water treatment with ozone | Optionally available |
| Dimensions [W x H x D] | 42 x 19 x 39 cm |
| Weight | 9.8 kg |
| Electrical data | 100 - 240 V AC; 50 - 60 Hz; 1.2 A; Protection class I |
| Classification | IIb medical device according to EU Directive 93/42/EEC |
| Warranty | 2 years in Europe |

| Form of therapy | OZOMED smartline |
|------------------------------------|------------------|
| Major autohemotherapy | ● |
| Minor autohemotherapy | ● |
| Rectal insufflation | ● |
| Bag treatment | ● |
| Low pressure ozone gas application | |
| Ozone injections | ● |
| Ozonized water | optionally |



Convenience kits for ozone therapy

OZONOSAN Set 4 Plus Convenience kit



- 1 OZONOSAN vacuum bottle 250 ml, microbubble system, 12 ml sodium citrate as anticoagulant, latex-free
- 1 Blood transfusion set with steel tip
- 1 Perfusion set optionally with 0.8 or 1.1 mm needle
- 1 Disposable syringe 50 ml, siliconized
- 1 Transfer filter set with cannula and filter (0.2 µm) for ozone introduction

OZONOSAN Set Rectal insufflation with dosing bellow



- 1 Dosing bellow
- 5 Ozone reservoirs with valve
- 50 Special disposable catheters 18 Ch.
- 50 Infusion lines 75 cm
- 50 Bacteria filter

Hänsler Medical GmbH

Nordring 8, 76473 Iffezheim, Germany

+49 (0) 7229 3046 0

@ info@haensler-medical.com

www.haensler-medical.com

