

OZONOSAN boardcase



Mobile ozone generator for home visits, as a second unit and much more.

The OZONOSAN boardcase is a fully mobile suitcase device with internal oxygen cylinder, e.g. for home visits. Excellent for all systemic and topical applications with ozone.



Product highlights

Continuous photometric ozone concentration measurement with automatic calibration.

Mobile unit, flexible and independent of patient location.

Operation with internal or externally connected oxygen cylinder.

Powerful catalyst for safe ozone destruction.

User-friendly operating concept.

Low maintenance with above-average product service life.

Developed and manufactured in Germany.

Technical specification	OZONOSAN boardcase
Device setup	Portable suitcase with wheels
Mode of operation	Normobaric infusion technique
Ozone concentration	1 - 50 µg/ml
Flow rate	0.6 - 10 l/min, depending on concentration
Operating pressure	0.5 bar
Vacuum pump	Optionally available, up to -0.8 bar adjustable
Dimensions [W x H x D]	55 x 19 x 40 cm
Weight	19.4 kg including oxygen cylinder
Electrical data	230 V AC; 50 Hz; 1.2 A; Protection class I
Classification	IIa medical device according to EU Directive 93/42/EEC
Warranty	2 years; 5 years for ozone tube and catalyst

Hänsler Medical - Only trust the original.

Form of therapy	OZONOSAN boardcase
Major autohemotherapy	●
Minor autohemotherapy	●
Rectal insufflation	●
Bag gassing	●
Low pressure ozone gas application	●
Ozone injections	●
Ozonized water	



Convenience kits for ozone therapy

OZONOSAN Set 4 Plus Convenience kit



- 1 OZONOSAN vacuum bottle 250 ml, microbubble system, 12 ml sodium citrate as anticoagulant, latex-free
- 1 Blood transfusion set with steel tip
- 1 Perfusion set optionally with 0.8 or 1.1 mm needle
- 1 Disposable syringe 50 ml, siliconized
- 1 Transfer filter set with cannula and filter (0.2 µm) for ozone introduction

OZONOSAN Set Rectal insufflation with dosing bellow



- 1 Dosing bellow
- 5 Ozone reservoirs with valve
- 50 Special disposable catheters 18 Ch.
- 50 Infusion lines 75 cm
- 50 Bacteria filters

Hänsler Medical GmbH

Nordring 8, 76473 Iffezheim, Germany

+49 (0) 7229 3046 0

@ info@haensler-medical.com

www.haensler-medical.com

