

The OZONOSAN photonik is ideal as a tabletop unit in smaller medical centers. This system generates ozone concentrations between 2 and $\leq 75 \mu\text{g/ml}$ for all common applications. A water ozonization unit is optionally available.



OZONOSAN photonik and Water Ozonization Unit

Product highlights

Continuous photometric ozone concentration measurement with automatic calibration.

Powerful catalyst for safe ozone destruction.

Water ozonation unit can be connected optionally.

User-friendly operating concept.

Low maintenance with above-average product service life.

Developed and manufactured in Germany.

Technical specification	OZONOSAN photonik
Device setup	Tabletop unit
Mode of operation	Normobaric infusion technique
Ozone concentration	2 - 75 $\mu\text{g/ml}$
Flow rate	0.8 - 8 l/min
Operating pressure	0.5 bar
Vacuum pump	Optionally available, up to -0.8 bar adjustable
Water treatment with ozone	OZONOSAN water ozonization unit optionally available
Dimensions [W x H x D]	20 x 56 x 37.5 cm
Weight	19,6 kg
Electrical data	230 V AC; 50 Hz; 1.2 A; Protection class I
Classification	Ila medical device according to EU Directive 93/42/EEC
Warranty	2 years; 5 years for ozone tube and catalyst

Form of therapy	OZONOSAN photonik
Major autohemotherapy	●
Minor autohemotherapy	●
Rectal insufflation	●
Bag gassing	●
Low pressure ozone gas application	●
Ozone injections	●
Ozonized water	optionally



Convenience kits for ozone therapy

OZONOSAN Set 4 Plus Convenience kit



- 1 OZONOSAN vacuum bottle 250 ml, microbubble system, 12 ml sodium citrate as anticoagulant, latex-free
- 1 Blood transfusion set with steel tip
- 1 Perfusion set optionally with 0.8 or 1.1 mm needle
- 1 Disposable syringe 50 ml, siliconized
- 1 Transfer filter set with cannula and filter (0.2 µm) for ozone introduction

OZONOSAN Set Rectal insufflation with dosing bellow



- 1 Dosing bellow
- 5 Ozone reservoirs with valve
- 50 Special disposable catheters 18 Ch.
- 50 Infusion lines 75 cm
- 50 Bacteria filters

Hänsler Medical GmbH

Nordring 8, 76473 Iffezheim, Germany

+49 (0) 7229 3046 0

@ info@haensler-medical.com

www.haensler-medical.com

